


# Two-Hole Osteotomy Support Plate System

Surgical Technique



## System Features

### Clinical Use


The 2-hole osteotomy support plate stabilizes hinge fractures that can occur unintentionally during an opening wedge osteotomy of the proximal tibia or distal femur.

### Implant System

The implant system contains a plate with two 4 mm × 30 mm nonlocking, self-tapping bone screws. The screw heads are nonthreaded and can be inserted in variable angles. The drill hole can be prepared with a 2 mm drill. To insert the screws, a 2.5 mm hexalobe screwdriver is needed.



## Surgical Technique Overview



**01**

### Identification of Fractured Bone


Perform a fluoroscopy in a strict A-P direction to identify the fractured cortical bone.



**02**

### Incision and Preparation of Bone


Incise the skin at the level of the fracture and open the deep fascia anterior to the iliotibial tract. Release the proximal attachment of the tibialis anterior muscle from the bone without stripping the periosteum. Retract the tibialis anterior muscle laterally to access the bone.



**03**

**Coverage of Fractured Bone With Placement of K-Wire**


Close the osteotomy gap to approximate the lateral bone and place the plate over the fracture. If desired, bend the plate to fit the anatomy. Secure the plate position with a 1.6 mm K-wire. Next, drill a 2 mm hole to a depth of 30 mm.



**04**

**Initial Screw Insertion**


Fully insert the first 4 mm self-tapping bone screw using a 2.5 mm hex screwdriver.



**05**

**Secondary Screw Insertion**

Remove the K-wire and drill a second hole of 30 mm depth, diverging away from the direction of the first screw. Insert the second bone screw.



**06**

**Final Implant Fixation**

Fully insert the second screw to complete the stabilization of the hinge fracture. The osteotomy can be reopened and continued to achieve the desired alignment correction.

## Ordering Information

### Instruments

2-Hole Osteotomy Support Plate Implant System	AR-13215
K-wire, 1.6 mm × 200 mm	AR-5050K-1
Drill bit, 2.0 mm, short, AO	AR-18800-18
Universal bending iron, osteotomy plates	AR-13322-02
Noncannulated short screwdriver shaft, 2.5 mm hex, Ø 5.5 mm × 134 mm long	AR-1995SHN

Products advertised in this surgical technique guide may not be available in all countries. For information on availability, please contact Arthrex Customer Service or your local Arthrex representative.

This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.



Arthrex manufacturer, authorized representative, and importer information (Arthrex eFUs)



US patent information