


Clavicle Plate and Screw System

Surgical Technique




Clavicle Plate and Screw System

The Clavicle Plate and Screw System is a comprehensive set of stainless steel plates and screws with instrumentation for treating central third and distal clavicle fractures. The plates are precontoured and anatomically designed to limit the need for bending, helping to ensure an optimal match with the patient's superior clavicle.



Achieve secondary fixation to the coracoid using the Knotless Distal Clavicle Plate Button TightRope® implant. The TightRope implant, consisting of a knotless clavicle plate button attached to a large pec button, is preloaded onto an inserter, allowing surgeons to perform an open, “push-through” technique without having to use a scope or access beneath the coracoid.



Clavicle Fracture Repair

Setup and Patient Positioning

Patient Positioning

Place the patient on the OR table. The beach chair position is recommended. Prep and drape the affected extremity in the normal sterile fashion. Place a roll or pad between the shoulder blades to allow retraction and aid in reduction. An arm holder, like the TRIMANO FORTIS arm holder, can also be very helpful in maintaining the position of the injured extremity.

Surgical Approach

Make a 3 cm to 5 cm horizontal incision over the superior clavicle. Subcutaneous dissection allows for identification of supraclavicular nerve branches.

Fracture Reduction


Reduce the fracture and use fluoroscopy to confirm reduction. Using reduction forceps can be very helpful in maintaining reduction. Apply interfragmentary (lag) screws as/if appropriate to help maintain fracture reduction. If necessary, smaller butterfly fragments may be reduced with 2.5 mm distal-tip fragment screws.

Plate Selection

Select the appropriate clavicle plate to match the patient's anatomy. Aluminum clavicle plate sizing templates may be used to help determine the appropriate implant. The plates are precontoured to reduce the need to bend. If contouring the plate is necessary, use the appropriate plate benders. Do not bend the plate near the locking holes as this may distort the threads and inhibit proper screw insertion. Repeated bending of the plate at the same location or by creating excessive acute angles may potentially lead to plate fatigue, failure, and/or breakage in situ.



TRIMANO FORTIS



Fragment Screws
2.5 mm Herbert screw design aids in reducing comminuted bone fragments



Clavicle Plate Sizing Template





Plate Placement

Place the plate onto the reduced clavicle and temporarily attach it to the bone using BB-Taks, K-wires, or plate-holding forceps.

Central Third Plates

| Product | Item Number |
|---|------------------------------|
|  | AR-2650CL Length = 86 mm |
|  | AR-2651CL Length = 93 mm |
|  | AR-2652CL Length = 98 mm |
|  | AR-2653CL Length = 120 mm |
|  | AR-2654CL Length = 87 mm |
|  | AR-2655CL Length = 76 mm |

| Product | Item Number |
|--|------------------------------|
|  | AR-2650CR Length = 86 mm |
|  | AR-2651CR Length = 93 mm |
|  | AR-2652CR Length = 98 mm |
|  | AR-2653CR Length = 120 mm |
|  | AR-2654CR Length = 87 mm |
|  | AR-2655CR Length = 76 mm |
|  | AR-2680ST Length = 76 mm |



- (a) Tapered end to minimize soft-tissue irritation
- (b) Locking holes
- (c) Compression slots that can accommodate plate buttons
- (d) Recessed slots for screws to sit flush with plate, minimizing soft-tissue irritation
- (e) Beveled edges to minimize soft-tissue irritation




Distal Clavicle Plates


| Product | Item Number |
|---|---------------------------------|
|  | AR-2622DL Length = 70 mm |
|  | AR-2656DL Length = 70 mm |
|  | AR-2657DL Length = 100 mm |
|  | AR-2685DL-10 Length = 126 mm |

| Product | Item Number |
|--|---------------------------------|
|  | AR-2622DR Length = 70 mm |
|  | AR-2656DR Length = 70 mm |
|  | AR-2657DR Length = 100 mm |
|  | AR-2685DR-10 Length = 126 mm |

| Product Description | Item Number |
|--|-------------|
| Knotless Distal Clavicle Plate Button TightRope® Implant | AR-2658TR |
| Compatible with any Arthrex distal clavicle plate stainless steel plate button; sits flush in any compression slot pec button; provides stable fixation beneath the coracoid | |



15° Divergent Screw Pattern




- (a) Tapered end to minimize soft-tissue irritation
- (b) K-wire hole for temporary plate fixation
- (c) Locking holes
- (d) Beveled edges to minimize soft-tissue irritation
- (e) Compression slots that can accommodate plate buttons


- (f) Threaded hole for positioning handle and distal drill guide
- (g) Chamfered suture eyelets to incorporate FiberWire® suture into the plate
- (h) 2.7 mm divergent locking holes
- (i) K-wire hole for temporary plate fixation

Central Third Plates


Nonlocking Screw Insertion



Place the 3.5 mm/2.5 mm drill guide into the appropriate plate slot and prepare a hole using the 2.5 mm drill bit. If drilling bicortically, place a retractor under the clavicle to help protect the neurovascular structures.




Measure for screw length using the screw depth device.




Select appropriate 3.5 mm or 4 mm screw and insert using the T15 hexalobe driver.


Locking Screw Insertion



Place the 3.5 mm threaded drill guide into the appropriate plate hole and prepare a hole using the 2.5 mm drill bit. If drilling bicortically, place a retractor under the clavicle to help protect the neurovascular structures. Read the corresponding screw length from the laser line on the drill. Alternatively, the screw depth device can be used to determine the screw length.




Select appropriate 3.5 mm locking screw and insert using the T15 hexalobe driver.



Insert remaining screws as needed to complete the repair. Arthrex BoneSync™ calcium phosphate cement or AlloSync™ Pure demineralized bone matrix can be used to help stabilize bone fragments within highly comminuted fractures and promote healing.

Distal Clavicle Plates


Locking Screw Insertion for 2.7 mm Screws




Place the drill guide onto the distal end of the plate so that the screw threads into the threaded plate hole and the pin seats into the K-wire hole.




Insert the guide sleeve into the desired hole. Leave the sleeve in place until after the screw is inserted.




Drill through the sleeve to the desired depth using the 2 mm drill. If desired, read the screw length on the drill where the line is level with the top of the guide sleeve.



Insert the screw depth indicator into the sleeve to determine the screw length.



Insert the 2.7 mm locking screw through the guide sleeve using the T10 hexalobe driver.




Alternate Method
Thread the 2.7 mm threaded drill guide into a 2.7 mm locking hole until fully seated. Using the 2 mm drill bit, drill to the desired depth and read the corresponding screw length from the laser line on the drill. Select the appropriate 2.7 mm screw and insert using the T10 hexalobe driver.


Coracoid Fixation

The clavicle and coracoid tunnels can be drilled either independently or fluoroscopy-assisted. A 10° to 15° cephalic Zanca view can optimize visualization of the coracoid and AC joint.


Plate the fracture using the appropriate distal clavicle plate. Identify the compression slot that allows the best anatomic placement for the Knotless Distal Clavicle Plate Button TightRope® implant. Do not place a screw in the designated compression slot.




1 Drill the clavicle tunnel using the 3.7 mm spade-tip drill. Fluoroscopy may be used to verify pin trajectory.




2 Drill the coracoid tunnel. Fluoroscopy may be used to confirm proper pin trajectory and drilling depth under the coracoid. A retractor can be placed under the coracoid to help prevent over-drilling the far cortex.



3 Insert the coracoid button into the clavicle tunnel. Light taps on the inserter handle can aid in advancing the button through the tunnels.




4 Insert the coracoid button through the clavicle and the coracoid tunnel. Fluoroscopy may be used to confirm proper button alignment and that it extends below the inferior cortex of the coracoid.




5

Remove the red pull tab from the handle to release the sutures and plate button from the inserter.




6

Turn the knurled stylet counterclockwise to release the coracoid button from the driver. With the inserter in the coracoid tunnel, grab the TightRope® implant sutures below the plate button and pull up to seat the button against the coracoid. Use fluoroscopy to verify that the button flipped properly. Remove the driver from the tunnels.



7

With the button firmly against the base of the coracoid, sequentially pull on the free suture limbs 1 cm to 2 cm at a time to reduce the button into the plate. A blunt hemostat may be placed under the plate button to aid in reduction.



8

Snap the button into the plate and use fluoroscopy to confirm reduction. Cut the free cinch limbs to complete the repair.

Confirm Reduction and Fixation

Confirm the final reduction and plate and screw fixation both visually and with fluoroscopy.

Plate and Screw Removal

If the plate and screws need to be removed, make an incision over the clavicle. Use the appropriate screwdriver to remove each screw. BoneSync™ calcium phosphate cement may be used to fill in the bone voids in the clavicle left from screw removal.

Post-op Protocol

Postoperatively, and until healing is complete, the fixation provided by this device should be considered temporary and it may not withstand weightbearing and/or other unsupported stress. To avoid adverse stress being applied to the device and help protect its fixation, the prescribed postoperative regimen should be followed strictly. Images should be obtained throughout the healing process to verify full bone healing before the patient is cleared for return to normal activities.

BONESYNC

Fast-Setting, Drillable Calcium Phosphate Cement


- BoneSync™ bone void filler is a fast-setting, drillable, resorbable cement composed of calcium phosphate and collagen.
- This synthetic bone graft sets within approximately 8 minutes and is a drillable solution for fracture repair.
- Once cured, BoneSync calcium phosphate cement can supplement hardware to support bone fixation and fill bony voids during surgical procedures.



ALLOSYNC PURE™

Demineralized Bone Matrix

- AlloSync™ Pure is a dehydrated osteoinductive demineralized bone matrix derived from 100% human allograft bone with no extrinsic carriers.
- AlloSync Pure bone matrix resists irrigation and can be used in a fluid environment (a).
- The proprietary rice-shape fiber technology used to process AlloSync Pure increases the osteoinduction and osteoconductive surface area to accelerate cellular ingrowth.



Ordering Information

Implants

| Product Description | Item Number |
|---|--------------|
| Clavicle Fracture Plate, central third, left | AR-2650CL |
| Clavicle Fracture Plate, central third, right | AR-2650CR |
| Clavicle Fracture Plate, central third, left | AR-2651CL |
| Clavicle Fracture Plate, central third, right | AR-2651CR |
| Clavicle Fracture Plate, central third, left | AR-2652CL |
| Clavicle Fracture Plate, central third, right | AR-2652CR |
| Clavicle Fracture Plate, central third, left | AR-2653CL |
| Clavicle Fracture Plate, central third, right | AR-2653CR |
| Fracture Plate, central third, left | AR-2654CL |
| Clavicle Fracture Plate, central third, right | AR-2654CR |
| Clavicle Fracture Plate, central third, left | AR-2655CL |
| Fracture Plate, central third, right | AR-2655CR |
| Clavicle Fracture Plate, 7 hole, straight | AR-2680ST |
| Distal Clavicle Fracture Plate, short, left | AR-2656DL |
| Distal Clavicle Fracture Plate, short, right | AR-2656DR |
| Clavicle Fracture Plate, long, left | AR-2657DL |
| Distal Clavicle Fracture Plate, long, right | AR-2657DR |
| Distal Clavicle Plate, 22°, left | AR-2622DL |
| Distal Clavicle Plate, 22°, right | AR-2622DR |
| Distal Clavicle Plate, 10 hole, left | AR-2685DL-10 |
| Distal Clavicle Plate, 10 hole, right | AR-2685DR-10 |

Coracoid Fixation

| Product Description | Item Number |
|---|-------------|
| Distal Clavicle Plate Button | AR-2658 |
| Acute AC Repair Kit | AR-2271 |
| Dog Bone™ Button | AR-2270 |
| Knotless Distal Clavicle Plate Button Assembly | AR-2658T |
| Knotless Distal Clavicle Plate Button TightRope | AR-2658TR |
| 3.7 mm Drill Pin | AR-2272 |

Clavicle Plate Set (AR-2650S)

| Product Description | Item Number |
|--|-------------|
| Cannulated Driver Handle | AR-13221AOC |
| BB-Tak | AR-13226 |
| Distal Clavicle Plate Positioning Handle | AR-2659 |
| Plate Bending Pliers | AR-2660 |
| Verbrugge Forceps w/ Pivoting Jaw | AR-2662 |
| K-wire, 0.062 in × 3 in | AR-2663 |
| T15 Hexalobe Driver | AR-8941DH |
| T10 Hexalobe Driver | AR-2665-T10 |
| Drill Guide, 2.6 mm/1.35 mm | AR-8943-03 |
| Bone Reduction Forceps | AR-8943-07 |
| Screwdriver, T10 Hexalobe | AR-8943-08 |
| Screwdriver, T15 Hexalobe | AR-8943-10 |
| Screw Holding Sleeve | AR-8943-11 |

Ordering Information

| Clavicle Plate Set (Cont.) | |
|--|-------------|
| Drill Bit, 2.5 mm | AR-8943-13 |
| Drill Guide, 3.5 mm/2.5 mm | AR-8943-14 |
| Depth Device | AR-8943-15 |
| Drill Bit, 2 mm | AR-8943-16 |
| Drill Guide, threaded, 2.7 mm | AR-8943-17 |
| Plate Bending Iron | AR-8943-18 |
| Freer Elevator | AR-8943-19 |
| Periosteal Elevator, 6 mm, curved blade | AR-8943-20 |
| Sharp Hook | AR-8943-21 |
| Hohmann Retractor, 15 mm | AR-8943-22 |
| Lobster Claw Forceps | AR-8943-23 |
| Drill Guide, threaded, 3.5 mm | AR-8943-26 |
| T15 Hexalobe Drive Shaft, long | AR-2666-T15 |
| Distal Clavicle Plate Drill Guide, left | AR-2664L |
| Distal Clavicle Plate Drill Guide, right | AR-2664R |
| Drill/Screw Guide Sleeve | AR-2668 |
| 2 mm Drill | AR-2669 |
| 2.7 mm Screw Depth Indicator | AR-2670 |
| 3.5 mm Locking Bending Guide | AR-8954-07 |
| Clavicle Plating System Instrument Case | AR-2650C |

| Clavicle Plate Sizing Templates | |
|-------------------------------------|--------------|
| Product Description | Item Number |
| Sizing Template for AR-2650CL/CR | AR-2650CT |
| Sizing Template for AR-2651CL/CR | AR-2651CT |
| Sizing Template for AR-2652CL/CR | AR-2652CT |
| Sizing Template for AR-2653CL/CR | AR-2653CT |
| Sizing Template for AR-2654CL/CR | AR-2654CT |
| Sizing Template for AR-2655CL/CR | AR-2655CT |
| Sizing Template for AR-2656DL/DR | AR-2656DT |
| Sizing Template for AR-2657DL/DR | AR-2657DT |
| Sizing Template for AR-2685DL/DR-10 | AR-2685DT-10 |
| Sizing Template for AR-2622DL/DR | AR-2622DT |
| Sizing Template for AR-2680ST | AR-2680STT |

| 2.7 mm Low-Profile Nonlocking Screws | |
|--------------------------------------|-------------|
| Product Description | Item Number |
| 2.7 mm Screw Caddy | AR-8827C-10 |
| Low-Profile Screw, 2.7 mm × 10 mm | AR-8827-10 |
| Low-Profile Screw, 2.7 mm × 12 mm | AR-8827-12 |
| Low-Profile Screw, 2.7 mm × 14 mm | AR-8827-14 |
| Low-Profile Screw, 2.7 mm × 16 mm | AR-8827-16 |
| Low-Profile Screw, 2.7 mm × 18 mm | AR-8827-18 |
| Low-Profile Screw, 2.7 mm × 20 mm | AR-8827-20 |
| Low-Profile Screw, 2.7 mm × 22 mm | AR-8827-22 |
| Low-Profile Screw, 2.7 mm × 24 mm | AR-8827-24 |

Ordering Information

2.7 mm Low-Profile Locking Screws

| Product Description | Item Number |
|---|-------------|
| Low-Profile Locking Screw, 2.7 mm × 8 mm | AR-8827L-08 |
| Low-Profile Locking Screw, 2.7 mm × 10 mm | AR-8827L-10 |
| Low-Profile Locking Screw, 2.7 mm × 12 mm | AR-8827L-12 |
| Low-Profile Locking Screw, 2.7 mm × 14 mm | AR-8827L-14 |
| Low-Profile Locking Screw, 2.7 mm × 16 mm | AR-8827L-16 |
| Low-Profile Locking Screw, 2.7 mm × 18 mm | AR-8827L-18 |
| Low-Profile Locking Screw, 2.7 mm × 20 mm | AR-8827L-20 |
| Low-Profile Locking Screw, 2.7 mm × 22 mm | AR-8827L-22 |
| Low-Profile Locking Screw, 2.7 mm × 24 mm | AR-8827L-24 |

3 mm Low-Profile Nonlocking Screws, cancellous

| Product Description | Item Number |
|--|-------------|
| Low-Profile Nonlocking Screw, 3 mm × 10 mm | AR-8830-10 |
| Low-Profile Nonlocking Screw, 3 mm × 12 mm | AR-8830-12 |
| Low-Profile Nonlocking Screw, 3 mm × 14 mm | AR-8830-14 |
| Low-Profile Nonlocking Screw, 3 mm × 16 mm | AR-8830-16 |
| Low-Profile Nonlocking Screw, 3 mm × 18 mm | AR-8830-18 |
| Low-Profile Nonlocking Screw, 3 mm × 20 mm | AR-8830-20 |
| Low-Profile Nonlocking Screw, 3 mm × 22 mm | AR-8830-22 |
| Low-Profile Nonlocking Screw, 3 mm × 24 mm | AR-8830-24 |

3.5 mm Low-Profile Nonlocking Screws, cortical

| Product Description | Item Number |
|---------------------------------------|-------------|
| Low-Profile Screw, SS, 3.5 mm × 8 mm | AR-8835-08 |
| Low-Profile Screw, SS, 3.5 mm × 10 mm | AR-8835-10 |
| Low-Profile Screw, SS, 3.5 mm × 12 mm | AR-8835-12 |
| Low-Profile Screw, SS, 3.5 mm × 14 mm | AR-8835-14 |
| Low-Profile Screw, SS, 3.5 mm × 16 mm | AR-8835-16 |
| Low-Profile Screw, SS, 3.5 mm × 18 mm | AR-8835-18 |
| Low-Profile Screw, SS, 3.5 mm × 20 mm | AR-8835-20 |
| Low-Profile Screw, SS, 3.5 mm × 22 mm | AR-8835-22 |
| Low-Profile Screw, SS, 3.5 mm × 24 mm | AR-8835-24 |
| Low-Profile Screw, SS, 3.5 mm × 26 mm | AR-8835-26 |
| Low-Profile Screw, SS, 3.5 mm × 28 mm | AR-8835-28 |
| Low-Profile Screw, SS, 3.5 mm × 30 mm | AR-8835-30 |

3.5 mm Low-Profile Locking Screws

| Product Description | Item Number |
|---|-------------|
| Low-Profile Locking Screw, 3.5 mm × 8 mm | AR-8835L-08 |
| Low-Profile Locking Screw, 3.5 mm × 10 mm | AR-8835L-10 |
| Low-Profile Locking Screw, 3.5 mm × 12 mm | AR-8835L-12 |
| Low-Profile Locking Screw, 3.5 mm × 14 mm | AR-8835L-14 |
| Low-Profile Locking Screw, 3.5 mm × 16 mm | AR-8835L-16 |
| Low-Profile Locking Screw, 3.5 mm × 18 mm | AR-8835L-18 |
| Low-Profile Locking Screw, 3.5 mm × 20 mm | AR-8835L-20 |
| Low-Profile Locking Screw, 3.5 mm × 22 mm | AR-8835L-22 |
| Low-Profile Locking Screw, 3.5 mm × 24 mm | AR-8835L-24 |

Ordering Information

4 mm Low-Profile Nonlocking Screws, cancellous

| Product Description | Item Number |
|--|-------------|
| Low-Profile Nonlocking Screw, 4 mm × 10 mm | AR-8840-10 |
| Low-Profile Nonlocking Screw, 4 mm × 12 mm | AR-8840-12 |
| Low-Profile Nonlocking Screw, 4 mm × 14 mm | AR-8840-14 |
| Low-Profile Nonlocking Screw, 4 mm × 16 mm | AR-8840-16 |
| Low-Profile Nonlocking Screw, 4 mm × 18 mm | AR-8840-18 |
| Low-Profile Nonlocking Screw, 4 mm × 20 mm | AR-8840-20 |
| Low-Profile Nonlocking Screw, 4 mm × 22 mm | AR-8840-22 |
| Low-Profile Nonlocking Screw, 4 mm × 24 mm | AR-8840-24 |

Fragment Screws

| Product Description | Item Number |
|--|-------------|
| Fragment Screw, 2.5 mm × 10 mm, hexalobe | AR-2665-10H |
| Fragment Screw, 2.5 mm × 12 mm, hexalobe | AR-2665-12H |
| Fragment Screw, 2.5 mm × 14 mm, hexalobe | AR-2665-14H |
| Fragment Screw, 2.5 mm × 16 mm, hexalobe | AR-2665-16H |
| Fragment Screw, 2.5 mm × 18 mm, hexalobe | AR-2665-18H |
| Fragment Screw, 2.5 mm × 20 mm, hexalobe | AR-2665-20H |
| Fragment Screw, 2.5 mm × 22 mm, hexalobe | AR-2665-22H |
| Fragment Screw, 2.5 mm × 24 mm, hexalobe | AR-2665-24H |

Required Fragment Screw Instruments

| Product Description | Item Number |
|----------------------|-------------|
| Step Drill | AR-2752 |
| T8 Hexalobe Driver | AR-2665-T8 |
| Small Holding Sleeve | AR-8943-45 |

TRIMANO FORTIS

| Product Description | Item Number |
|--------------------------------------|-------------|
| TRIMANO FORTIS Support Arm | AR-1740 |
| TRIMANO FORTIS Adapter | AR-1741 |
| TRIMANO®* Arm Holder Beach Chair Kit | AR-1644 |

Orthobiologics

| Product Description | Item Number |
|----------------------|-------------|
| BoneSync™ Kit, 5 cc | ABS-3105 |
| AlloSync™ Pure, 5 cc | ABS-2010-05 |

*TRIMANO is a registered trademark of Maquet GmbH.



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.

View U.S. patent information at www.arthrex.com/corporate/virtual-patent-marking

© 2019 Arthrex, Inc. All rights reserved. | www.arthrex.com | LT1-0255-EN_C