


Collagen Coated FiberWire[®] Suture

FiberWire is now available with a Type-I bovine collagen coating to offer our customers a softer coating for improved handling characteristics.




130x



Collagen
Coated FiberWire



SpeedFix[™] with Collagen Coated FiberTape




SpeedBridge[™] with Collagen Coated FiberTape




SpeedWhip[™]
with Collagen Coated FiberWire




Collagen Coated FiberTape® SpeedBridge Repair




A Collagen Coated FiberTape and a standard FiberTape were used in a cadaveric SpeedBridge repair. This test proved that the collagen remains intact throughout a procedure. Both of these FiberTapes were dipped in a collagen dye prior to testing.



Enhanced zoom 130x



Standard #2 FiberWire



Collagen Coated #2 FiberWire

Ordering Information:

Collagen Coated #2 FiberWire w/ tapered needle
Collagen Coated #2 FiberLoop w/ straight needle
Collagen Coated FiberTape, 7"

AR-7200B
AR-7234B
AR-7237-7B



Scan to view
more product
information

<http://tinyurl.com/CCFIBERWIRE-2013>



Arthrex®

©2013, Arthrex, Inc. All rights reserved.
LS1-0235-EN_A